

# EZZI-LIFT®

## PRODUCT SPECIFICATIONS



The Ezzi-lift microcurrent device for aesthetic use gently stimulates, tightens, and tones the face and body for a more youthful, natural appearance. The non-invasive, hand-held device can be used without creams or gels and is ideal for use at home or in the clinic.

### BENEFITS

- Diminishes the appearance of fine lines and wrinkles
- Tightens skin around the neck and jawline
- Delivers gentle microcurrent massage to the face and body
- Enables treatment without the need for creams or gels

### FEATURES

- Avazzia BEST® microcurrent platform
- Medical-grade stainless steel electrodes
- Multiple treatment modes and intensity options
- Accessory port for optional treatment tools

## TECHNICAL SPECIFICATIONS

Size, weight	<b>Size:</b> 119mm X 64mm X 33mm (4.69 in. x 2.52 in. x 1.30 in.)   <b>Weight (device only):</b> Approx. 108 grams (3.8 oz)   <b>Weight (with batteries):</b> Approx. 153 grams (5.4 oz)
Power source	2 AA 1.5V alkaline batteries
Pulses, output	<b>Pulse Duration:</b> less than 2ms   <b>Frequency of outputs:</b> 15 to 121 Hz <b>Maximum Voltage:</b> less than 42V at 500Ω   <b>Maximum Current:</b> less than 500μA
Waveform	Pulsed current, damped, biphasic, asymmetric sinusoidal with microcurrent output

## ADDITIONAL PRODUCT INFORMATION

Modes	Face Stimulate, Face Smooth, Eye Stimulate, Eye Smooth, 60-minute Timeout
User Controls	ON/OFF, Power output intensity increase (▶) and decrease (◀), Mode selection
Indicators: audible, visual	<b>Mode selection:</b> LED indicates mode selection. <b>Power intensity:</b> LEDs indicate power level. <b>Time:</b> Chirps and beeps. A quiet, subtle tick sound occurs every 30 seconds. A double beep occurs every two minutes. <b>Low battery condition:</b> Indicated by a long series of twenty beeps, in concert with blinking mode LEDs, and followed by device shut off.
Available accessories	Y-electrode, brush electrode, pencil electrode, lead wires (SOLD SEPARATELY)
Availability, clearances	US FDA cleared aesthetic use over-the-counter, CE approved, CB certified, Health Canada certified Cosmetic device available for general purchase. Visit <a href="http://www.avazzia.com">www.avazzia.com</a> for information.
Quality manufacturing	Assembled by Avazzia in U.S. with imported and domestic components. ISO 13485 international certification

WARNING: NOT TO BE USED BY THOSE WITH PACEMAKER OR OTHER IMPLANTABLE ELECTRICAL DEVICE. DO NOT USE DIRECTLY OVER CAROTID ARTERY, DIRECTLY ON OPEN WOUNDS, OR CAUSE CURRENT TO RUN TRANSCRANIALY. SEE OWNER'S MANUAL FOR COMPLETE SAFETY AND USE INSTRUCTIONS.