

AVAZZIA BLUE™ DEVICE

Pain Relief Without Drugs

Avazzia Blue technology is the same advanced neuromodulator that is built into prescription devices for pain relief including the same unique complex waveform with pulsed, high-voltage microcurrent pulse streams.

Interactive chirp and ring indicate the tissue is changing or if the tissue has stopped changing.

The US FDA has reviewed and cleared the Avazzia Blue™ OTC device for the temporary relief of pain associated with sore and aching muscles in the shoulder, waist, back, back of the neck, upper extremities (arm), and lower extremities (leg) due to strain from exercise or normal household work activities.

"After spending eight hours a day using a computer mouse, I began suffering from severe numbness in my right hand and arm. I was introduced to the Avazzia device and began seeing improvement almost immediately. Within six months I was totally symptom-free and able to avoid surgery (which I was told was my other option)."

--- Dianne Stultz, Plano, TX



Non-Invasive &
Non-Drug Pain Relief

Complex Microcurrent
Electrotherapy
Pulse Streams

Automatic Interactive
Cybernetic Bio-Electric
Feedback

Advanced Neurological
Signals to Maximize
Pain Relief

Handheld &
Battery Operated

Simple and Easy To Use

No Prescription Required

Providing pain relief without drugs or surgery since 2004 in Dallas, Texas, Avazzia designs, manufactures, markets and distributes scientifically advanced, medical devices for management of pain.

Avazzia products incorporate patented B.E.S.T.™ technology which is Bio-Electric Simulation Technology. BEST technology senses the body's responses to the neurostimulation and detects, measures, analyzes and responds to the body.

- ISO 13485 certified
- FDA cleared for OTC pain relief
- HealthCanada
- CE certified
- CB certified

AVAZZIA™
Innovation in Health & Wellness Technology

13140 Coit Road, Suite 515
Dallas, Texas 75240
tel 214.575.2820
www.avazzia.com

AVAZZIA™
Innovation in Health & Wellness Technology

AVAZZIA BLUE™ OTC DEVICE

Prescription level pain relief without drugs



Easy To Use

The Avazzia Blue™ OTC BEST medical device incorporates two operating modes per device accessed by simply pressing the MODE switch.



Power levels are accessed by pressing “+” and “-” switches. Simple, easy to see LEDs that indicate different modes and low battery indicator along with automatic shutoff are standard features.



NOTE:
BEST™ devices are even more convenient with compatible conductive garments.

Directions for Use

1. Turn on the device.
2. Connect the device to the affected area by either:



- a. placing the on-board electrodes firmly on skin

or



- b. attaching electrodes to skin and then to device.

3. Select a comfortable power level.
4. Select an operating mode.
5. Apply for several minutes.
6. Select the next mode and repeat.

Complete Solution to knock out pain

ACCESSORIES

Sold Separately

Self-adhesive electrode pads
Standard with kit.



Y Electrode



For hard to reach areas

Videos Online

Free device instructional videos are available online at:
AvazziaTraining.com under the products tab for Avazzia Blue.

AVAZZIA
Innovation in Health & Wellness Technology

- Advanced neuromodulation signals maximize pain relief.
- Can be used by itself or with accessories.
- 2 AA battery operated.
- Manufactured in Texas and sold globally
- Over-the-counter
- FDA 510k cleared microcurrent TENS for OTC pain relief.