


Date: \_\_\_\_\_

Patient ID Number: \_\_\_\_\_

Site Location: \_\_\_\_\_

Physician/Investigator: \_\_\_\_\_



Treatment notes:

- Follow Little Wings Treatment directions.
- Initial Reading may be taken at either the sternocleidomastoid muscle (SCM) on the side of the neck or the trapezius muscle at the top of the shoulder. In all cases, avoid the carotid artery and jugular vein.
- Take readings in the same location (SCM or trapezius muscle) on both sides of the body.
- IR reading is taken in Relax/Assess mode (default mode); switch to Acute mode for Little Wings therapy. Power setting and duration can be read from the PRO-SPORT screen after Little Wings treatment is complete. Record those numbers. Then switch back to Relax/Assess for final IR.
- Readings and therapy may be done with onboard electrode or with Y-electrode attachment.

	Pre-therapy IR reading (Relax/Assess)	Power setting in Acute	Duration to achieve five contractions	IR reading (Relax/Assess) after Little Wings
Left Side				
Right Side				

**Treatment Notes**