



## BEST-PRO 1 Product Description

The BEST™ line of products are hand-held, battery-operated, electro-stimulation medical devices. The Avazzia BEST-PRO 1 is designed for at-home use and offers easy-to-use features including four popular frequencies that allow patients to administer microcurrent pain relief at home.

Research has shown microcurrent offers long-term pain relief without drugs or surgery. It is not habit-forming and non-addictive.



PRODUCT SPECIFICATIONS	
SIZE AND WEIGHT	<b>Size:</b> 117mm X 64mm X 33mm. <b>Weight (device only):</b> Approx. 128 grams (4.5 oz). <b>Weight (with batteries):</b> Approx. 173grams (6.1 oz)
FREQUENCY, HZ	4 Hz to 352 Hz
POWER SOURCE	2 AA 1.5V DC batteries
OUTPUT	<b>Output Voltage:</b> 20 V to 500 V <b>Output Current:</b> 2 mA to 90 mA <b>Current Density:</b> 0.1 to 966 $\mu\text{A}/\text{cm}^2$ <b>Power Density:</b> 0.1 to 1.2 $\text{mW}/\text{cm}^2$
PULSES	<b>Pulse Duration:</b> less than 2 ms
WAVEFORM	Pulsed, damped, biphasic sinusoidal
AUTOMATIC SHUT-OFF	60 minutes
PRODUCT INFORMATION	
MODES OF OPERATION	Assess Relax; Stimulate; Deep Stimulate; Acute
REACTIONS	In Assess Relax Mode
AUDIBLE INDICATIONS	30 second tick, 2 minute double beep, chirp in Assess Relax mode at start of application, ring in Assess Relax mode indicating progress
VISUAL INDICATIONS	Mode selection and power intensity indicating LEDs
PRODUCTION	Manufactured by Avazzia in U.S.
AVAILABILITY	By prescription only in the U.S. For assistance in filling prescription, contact Avazzia.
ACCESSORIES (SOLD SEPARATELY)	Y-electrode, adhesive conductive pads, brush electrode, pencil electrode
CLEARANCES	ISO 13485 international certification CE approved, CB certified HealthCanada certified FDA-Cleared for Pain
WARNINGS	Not to be used by those with pacemaker or other implantable electrical device. Not to be used by pregnant or nursing women. Do not use directly over carotid artery or cause current to run transcranially. Do not use directly on open wounds.