



Avazzia Blue™ Product Performance Specifications

This product is a portable, hand-held, battery-operated, microcurrent electro-stimulation medical device. The Avazzia Blue offers an easy-to-use microcurrent device featuring two of the effective frequencies, Blue Relax and Blue Stimulate. The simplicity of design makes it a perfect choice for patient use at home.

Avazzia BEST technology uses unique output waveforms and frequencies and treatment will not be accommodated by the body. The microcurrent treatment may be administered using the onboard electrodes on the unit's back or with accessories (purchased separately).



	PRODUCT PERFORMANCE SPECIFICATIONS
DEVICE SPECIFICATIONS	Size: 117mm X 64mm X 31mm. Weight (device only): Approx. 108 grams (3.8 oz). Weight (with batteries): Approx. 153 grams (5.4 oz)
FREQUENCY, HZ	25 to 200 Hz
POWER	Power Source: 2 AA 1.5V DC batteries
OUTPUT	Maximum Voltage: up to 650V Maximum Current: 1300µa
PULSES	Pulse Duration: less than 2 ms
WAVEFORM	Pulsed, damped, biphasic sinusoidal
MODES	Blue Relax (indicator light blinks slowly) Blue Stimulate (indicator light blinks rapidly)
REACTIONS	D and ZERO Reactions in Blue Relax mode (indicated by audio and visual cues).
AUTOMATIC SHUT-OFF	60 minutes

	PRODUCT INFORMATION
MODES OF OPERATION	Blue Relax and Blue Stimulate
CLEARANCES	ISO 13485 international certification CE approved, CB certified HealthCanada certified FDA 510k Cleared for OTC relief of pain associated with sore and aching muscles in the shoulder, waist, back, back of the neck, upper extremities (arm), and lower extremities (leg) due to strain from exercise or normal household work activities.
PRODUCTION	Assembled by Avazzia in U.S. with imported and domestic components
AVAILABILITY	Medical device available for general purchase. Visit www.avazzia.com for information.
ACCESSORIES (SOLD SEPARATELY)	Y-electrode, adhesive conductive pads, lead wires
CONTRAINDICATIONS	Not to be used by those with pacemaker or other implantable electrical device or who may be pregnant or nursing . See owner's manual for safety information.